
Data Sheet

GFP Selector

Catalog No.	N0310-S	N0310-M	N0310-L
Resin μ L	100	250	1000
Slurry μ L	200	500	2000

Description

GFP Selector is based on a high-affinity single-domain antibody (sdAb) that is covalently immobilized on 4 % cross-linked agarose beads. The innovative, oriented and selective attachment via a flexible linker guarantees a high accessibility of the sdAbs and largely eliminates batch-to-batch variations. Due to the single-chain nature of sdAbs and their stable and covalent attachment, no leakage of light and heavy chains is observed during elution with SDS sample buffer. GFP Selector thus features high affinity and superior capacity for GFP fusion proteins while showing negligible unspecific background. GFP Selector immobilizes a wide range of GFP derivatives.

GFP Selector is compatible not only with physiological buffers but also with high stringency buffers (see "Buffer Compatibility" below). GFP Selector thus provides great freedom to adjust the binding and washing conditions to the experimental needs.

For recommended protocols please see our webpage at www.nano-tag.com/protocols.

**To be used *in vitro* / for research only,
not for diagnostic or therapeutic use!**

Non-toxic, non-hazardous, non-infectious.

Properties

Support 4 % cross-linked agarose, bead size 50-150 µm

Coating sdAb anti-GFP clone 1H1

Reactivity Recognizes GFP, mEGFP, superfolder GFP and most common CFP and YFP variants.

Does not cross-react with mCherry, mRFP, dsRed, mTagBFP, mTagRFP or their most common derivatives.

For a detailed overview visit our webpage at www.nano-tag.com/protocols.

Capacity > 4 µg GFP per µl of packed beads

Formulation 50 % slurry in PBS containing 20 % Ethanol

Shipment Shipped at ambient temperature

Storage Store at 4 °C, do not freeze

Stability Stable for 12 months

Buffer Compatibility

- Common buffer substances at pH 5 to 9
- 2% Triton X-100, 1% Tween-20, 1% NP-40, 1% CHAPS, 1% Deoxycholate, 0.1% SDS
- 4 M NaCl, 2 M KCl, 1 M MgCl₂, 100 mM EDTA
- 4 M urea
- 10 mM DTT, 10 mM 2-Mercaptoethanol
- RNase A, DNase I, Benzonase, protease inhibitors

Limited Use Label License

The purchaser of this product is granted a limited, non-transferable right to use the purchased material only for internal, academic research purposes. The license explicitly excludes the right to resell this product or any of its components. For information on further licensing possibilities please contact info@nano-tag.com

Fluorescent Protein	GFP Selector	RFP Selector	TagFP Selector
Sirius	✓		
tSapphire	✓✓		
Cerulean	✓✓✓		
eCFP	✓✓✓✓		
mTurquoise, mTurquoise2	✓✓✓✓✓		
acGFP	✓✓✓✓✓		
EGFP	✓✓✓✓✓		
Emerald GFP	✓✓✓✓✓		
mEGFP (A206K)	✓✓✓✓✓		
mEGFP (L221K)	✓✓✓✓✓		
superecliptic pHluorin	✓✓✓✓✓		
paGFP	✓✓✓✓✓		
superfolder GFP	✓✓✓✓✓		
eYFP	✓✓✓✓✓		
mVenus	✓✓✓✓✓		
Citrine	✓✓✓✓✓		
mRFP		✓	
mCherry		✓✓	
mScarlet		✓✓✓	
mOrange2		✓✓✓✓	
tdTomato		✓✓✓✓✓	
dsRed1 , dsRed2		✓✓✓✓✓	
mTagBFP			✓
mKate, mKate2			✓✓
mTagRFP			✓✓
mTagRFP657			✓✓

Not recognized: Dendra2, Dronpa, tdEOS, mEOS3.2, mRuby3 and mTFP (mTeal)