

Data Sheet

FluoTag-X2 anti-Chicken IgY Fluorophore-conjugated

Fluorophore	Catalog No. (500µl)	Catalog No. (2000µl)
Abberior-STAR RED	N0702-AbRED-S	N0702-AbRED-L
Abberior-STAR 580	N0702-Ab580-S	N0702-Ab580-L
Abberior-STAR 635P	N0702-Ab635P-S	N0702-Ab635P-L
AlexaFluor 647	N0702-AF647-S	N0702-AF647-L
Atto488	N0702-At488-S	N0702-At488-L
Atto542	N0702-At542-S	N0702-At542-L
Atto565	N0702-At565-S	N0702-At565-L
Atto647N	N0702-At647N-S	N0702-At647N-L
Sulfo-Cyanine 3	N0702-SC3-S	N0702-SC3-L
Sulfo-Cyanine 5	N0702-SC5-S	N0702-SC5-L

Identity Camelid sdAb anti-Chicken IgY

Content 5 µM fluorescently labeled sdAb in buffered saline, 50% glycerol, 0.09% sodium azide.

Storage Up to 3 months: -20 °C
Up to 12 months: -80 °C or below
Protect from light!

Applications IF / ICC / IHC Recommended dilution 1:500

Specificity Recognizes specifically Chicken IgY
Does not cross-react with antibodies from species like: Mouse, Rabbit, Guinea Pig, Rat or Goat.

For more information please visit our web page at www.nano-tag.com

To be used *in vitro* / for research only.
Non-toxic, non-hazardous, non-infectious.