

## Data Sheet

### FluoTag-X2 anti-Mouse IgG1 Fluorophore-conjugated

Fluorophore	Catalog No. (500 µl)	Catalog No. (2000µl)
Abberior-Star 580	N2002-Ab580-S	N2002-Ab580-L
Abberior-Star RED	N2002-AbRED-S	N2002-AbRED-L
Abberior-Star 635P	N2002-Ab635P-S	N2002-Ab635P-L
AlexaFluor647	N2002-AF647-S	N2002-AF647-L
Atto488	N2002-At488-S	N2002-At488-L
Atto542	N2002-At542-S	N2002-At542-L
Atto565	N2002-At565-S	N2002-At565-L
Atto647N	N2002-At647N-S	N2002-At647N-L
Sulfo-Cyanine 3	N2002-SC3-S	N2002-SC3-L
Sulfo-Cyanine 5	N2002-SC5-S	N2002-SC5-L
Licor 680RD	N2002-Li680-S	N2002-Li680-L
Licor 800CW	N2002-Li800-S	N2002-Li800-L

**Identity** Camelid sdAb anti-Mouse IgG1, Clone 10A4

**Content** 5 µM fluorescently labeled sdAb in buffered saline, 50% glycerol, 0.09% sodium azide.

**Storage** Up to 3 months: -20 °C  
Up to 12 months: -80 °C or below  
Protect from light!

**Applications** IF/ICC/WB Recommended dilution 1:500

**Specificity** Specifically recognizes Mouse IgG1 isotype.  
Does not cross-react with antibodies from species like: Rat, Guinea Pig, Rabbit, Chicken or Goat.

For more information please visit our web page at [www.nano-tag.com](http://www.nano-tag.com)

To be used *in vitro* / for research only.  
Non-toxic, non-hazardous, non-infectious.